


## GETTING STARTED WITH THE INTERX DEVICE

We may call this a “device,” but the InterX is 510(k) cleared by the U.S. FDA for pain relief and pain management, carries the European CE mark, and will serve as an important tool in your pain management toolbox. To help familiarize you with the InterX and its capabilities, here are a few things you should know.



The InterX 1000's built-in timer clicks and flashes a light every 30 seconds to help you time treatments. So, for example, if the protocol calls for 2 minutes of treatment, hold the device on the area for 4 clicks. The device is designed to shut off after 2 minutes of no use. See the device manual for more instructions and use.

**To preserve the life of the device, only use high-quality batteries.**





- Preset 1** = 15-60hz
- Preset 2** = 60hz
- Preset 3** = 30-120hz
- Preset 4** = 240hz
- Preset 5** = 480hz

Use a lower setting for patients new to the InterX 1000. On the first day of a new treatment program, start on Preset 1, the lowest frequency setting. Use this setting for the first 1-2 days then switch to Preset 2 for the next 1-2 days of treatment.

As the body gets used to the stimulation, go to Preset 3, then Preset 4. You'll want to rotate through Presets settings 1-4 throughout the treatment program.

It helps to make a note of which Preset is used every day and what areas are treated.

Do not progress to the next highest preset setting if symptoms increase at a lower preset setting. Keep treating at the lower settings until positive results are seen.

Preset 5 is for acute or emergency treatments or very strong pain that needs to be worked through. Please note, setting 5 can either work very well on strong pain or irritate it. By design, setting 5 does not work on the flex arrays.

INTERX

**If you need additional information on the device or its use, please refer to the supplied Operations Manual or call our Support Team at (972) 665-1810.**