



SCALP TREATMENT PROTOCOL

Before beginning this treatment protocol, please read the “Getting Started” documents to familiarize yourself with the InterX 1000 and the treatment methods.


Scan and Target



Use the comb attachment electrode for this application.

Turn the device on and select the desired preset and intensity. For the first 1-2 treatments, use setting one or two.

Place the electrode onto the area around the anterior hairline to set the intensity. This ensures that you have good skin contact for the initial sensation to be correct.



1.

Scalp treatment helps reduce neuropathic pain. If the pain originates on the right side of the body, for example, in the right hand, you can treat the left side of the head in the scalp. In other words, the right side of the brain controls the left side of the body and vice versa.



Comb the electrode very slowly through your hair two or three times from the anterior hairline to the posterior hairline.

2.

Areas in the scalp that are sharp or more tender must be treated more. This is the area of the brain most affected. Some of the worst pain is reduced by treatment in the scalp using the comb electrode. Suggested time of treatment is 10-15 minutes.

Target


Target those points where you feel stronger stimulation. Target any sensitive points for 1 minute each. You may need to reduce the intensity while treating these points.

Continue like this until the whole scalp has been scanned and active points targeted.

You may have a delayed sensation, in which case it will be necessary to hold the comb in one spot for 10 seconds before the stimulation penetrates the scalp. This delayed response will shorten with subsequent treatments.

3.

Make sure you maintain good contact between the comb electrode and the skin on the scalp when treating over the scalp, covering the whole scalp in columns and moving slowly.



Anterior hairline

Posterior Hairline