




FACE TREATMENT PROTOCOL

Before beginning this treatment protocol, please read the “Getting Started” document and familiarize yourself with the InterX 5002. If you have any questions about these instructions, some may be answered in the operations manual that comes with the device, or contact your sales representative.

When there is chronic pain, the human body often remains in a sympathetic state, better known as a fight or flight response, which dictates how the body perceives and reacts to a threat. But the body has difficulty shifting to the parasympathetic state, a state where the body can rest, relax, digest, and heal.



Treating the face is a relaxing and necessary treatment that produces results. It can serve as the catalyst for the healing process by helping to move the body into a parasympathetic state.

1.

When the device is scanned or applied over a strong or active point, it might present as sticky, sharp, or quiet.

2.

To maximize results, increase treatment time to 4 minutes per point, for a total of 32 minutes, scanning the electrode toward the ear at each level.

3.


If the patient experiences sleeplessness, anxiety, or PTSD, treat the face before any other treatment.

4.

Remember to focus extra time on the level of the face related to the patient's pain.

Treating the face benefits the entire body

The areas in figure 1 marked V1, V2, and V3 are known as dermatome areas and can affect symptoms throughout the body.



Treat the face at least once a day

Lightly dampen the skin to facilitate stimulation.

Begin by selecting 15hz. After treating the face several times, switch to 60, 180, 240, then rotate the settings on days of treatment. Settings 15 and 60 are calming; 180 and 240 may be energizing.

Apply the built-in electrode or dome to each of the six points indicated on the face in figure 1. Make sure the surface of the electrode provides complete contact with the skin.

Standard treatment time for each point on the face is two minutes or four clicks of the device's internal timer. You can treat longer for a greater effect. See Pro Tip #2.